

Interval XXL

Care bed - 120 / 140 cm

The 3-functions Interval XXL bed has electrical adjustments and is available in 120 and 140 cm. Delivered in a white or beige frame with an extra strong, grid mattress base. The legpart comes with an electrically adjustable knee-buckle whereas the footpart can be adjusted manually. The complete bed can be ordered with wooden or metal side rails.



Comfortable and practical

The Interval XXL bed is equipped with an electrically adjustable head part and knee buckle. With the simple touch of a remote button, you can position the bed to read, watching TV or dine. The leg rest in its turn can also be adjusted to the most comfortable position.

This variable height bed is equipped with 4 wheels with separate brakes, ensuring its stability at all times.

Aesthetic and easy to maintain

The foot and head board are an integral part of this XXL bed. The warm colour of the wood makes this bed an elegant piece of furniture that will fit in almost any interior.

Frame colors



Easy-to-use remote with safety locks



Metal mattress holders on all sides



Manual adjustment of the leg part



Easy handling knobs to lower or raise the side rails

Technical data (120 cm mattress bed)

| | |
|---|-----------------------------|
| ■ Total width / length / height | 1330 x 2080 x 860 - 1295 mm |
| ■ Total height of mattress support | 370 - 805 mm |
| ■ Mattress spec. | 1200 x 2000 x 120 mm |
| ■ Weight | 124 kg |
| ■ Max. load (incl. patient, mattress and options) | 275 kg |
| ■ Additional info | CE |

Technical data (140 cm mattress bed)

| | |
|---|-----------------------------|
| ■ Total width / length / height | 1530 x 2080 x 860 - 1295 mm |
| ■ Total height of mattress support | 370 - 805 mm |
| ■ Mattress spec. | 1400 x 2000 x 140 mm |
| ■ Weight | 139 kg |
| ■ Max. load (incl. patient, mattress and options) | 275 kg |
| ■ Additional info | CE |

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer