

Four-Light

R O L L A T O R



Width	660 mm
Length (extended)	750 mm
Seat width	450 mm
Seat height	550 mm
Brakes	multifunctional brake
Max. inclination	(forward) 15° (rear and sideways) 10°
Weight	8 kg
Weight with basket	8,45 kg
Max. users weight	130 kg
Max. content basket	5 kg
Wheels	front: Ø 200 mm x 35 mm rear: Ø 175 mm x 35 mm
Colours	carbon grey - orange - white

Deviations in sizes and weight
(15 mm and 1,5 kg) are possible



N.V. Vermeiren N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: +32 (0)3 620 20 20 - Fax: +32 (0)3 666 48 94 - info@vermeiren.be - website: www.vermeiren.com