



## BATTERY DECLARATION

**The manufacturer:**

VERMEIREN GROUP

Actor ID/SRN: BE-MF-000019750

**Address:**

Vermeirenplein 1/15

2920 Kalmthout

Belgium

**declares under his sole responsibility for the CE marked medical devices:**

Product group: Electric (Powered) wheelchairs, Scooters

Brand: VERMEIREN

That the batteries used in our electric (powered) wheelchairs, and scooters are 2 x 12V lead acid batteries with a capacity ranging from 12Ah to 150Ah (288Wh to 3600Wh) depending on type and model from the following manufacturers:

- Cell Power
- Dyno Europe
- MK
- E-drive
- Long
- Celectric
- Tianneng

They are marked as NON-SPILLABLE (UN 2800) and they are EXEMPT from the international hazardous materials (also known as "dangerous goods") regulations since they comply with the following requirements:

- The vibration and pressure differential tests found in Packing Instruction 806 and Special Provision A67 of the **International Air Transport Association (IATA) Dangerous Goods Regulations**;
- The vibration and pressure differential tests found in Packing instruction 806 and Special Provision A67 of the **International Civil Aviation Organization (ICAO) Technical Instructions for the Safe Transport of Dangerous Goods by Air**;
- The vibration, pressure differential, and "crack" tests found in Special Provision 238 of the **International Maritime Dangerous Goods (IMDG) Code**.
- The vibration, pressure differential, and "crack" tests found in Special Provision 238 and 598 of the **International Carriage of Dangerous Goods by Road (ADR) Code**.

and therefore are unrestricted for transportation by any means.

City, Date: Kalmthout, 2025-04-09

Signature:

Name: Patrick Vermeiren

Function: CEO

**VERMEIREN GROUP**

Vermeirenplein 1-15

2920 Kalmthout – Belgium

Tel. +32 3 620 20 20

VAT BE 0887.199.216 - RPR Antwerpen

[www.vermeiren.com](http://www.vermeiren.com) | [info@vermeiren.com](mailto:info@vermeiren.com)

**VERMEIREN**  
we care for you

Certified ISO 13485