



DECLARATION ON ISOLATION OF ELECTRICAL CIRCUITS

The manufacturer:
VERMEIREN GROUP
Actor ID/SRN: BE-MF-000019750

Address:
Vermeirenplein 1/15
2920 Kalmlthout
Belgium

declares under his sole responsibility for the CE marked medical devices:

Product group: Electric (Powered) wheelchairs, Scooters
Brand: VERMEIREN

That, for those products indicated as 'suitable for transport by aircraft', electrical circuits can be isolated according to one of the following methods (in order of preference) as proposed by the "IATA Guidance Document for Battery-Powered Wheelchairs and Mobility Aids – 2026" :

- 1) Lock the battery in place, and switch off the battery (either on the battery itself, or with the battery switch where available);
- 2) Disconnect/Unplug the joystick module;
- 3) "Lock out" the medical device by activating the console lock.

Refer to the instructions for use to determine the most suitable method for your product.

City, Date: Kalmlthout, 2026-09-07
Signature:

Name:
Function:


VERMEIREN Group NV
ing. Tom Stal, MSc
Compliance & Regulatory Affairs
Vermeirenplein 1-15
2920 Kalmlthout, Belgium
QA/RA

VERMEIREN GROUP
Vermeirenplein 1-15
2920 Kalmlthout – Belgium
Tel. +32 3 620 20 20
VAT BE 0887.199.216 - RPR Antwerpen
www.vermeiren.com | info@vermeiren.com

VERMEIREN
we care for you

Certified ISO 13485