


# Express

## ELECTRONIC FOLDABLE WHEELCHAIR

Electronic wheelchair with a double cross, equipped with footrests and armrests adjustable in height and width. Equipped with a powerful engine of 2 x 150W and 38 A/h batteries.



Easy and user-friendly controller

## FEATURES

Disassemble fast and easy

### LINX™ ACCESS

#### Wheelchair Control System

LE Controller System :

- Ideal for light weight, indoor wheelchairs
- Easy to use
- Simple
- Reliable



Very compact and foldable chair which facilitates transport



2 separate batteries

Wide range of accessories and options.



2x12V-38Ah 2x12V-50Ah	2x150W	6 km/h	18 km 25 km	7°	15	1400
--------------------------	--------	--------	----------------	----	----	------



580		390								
610		420								
630		440								
650	400	460	225 - 315	540	430	430	400-520	950	67	90
690		500								

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

