

V200 D

STANDARD & LIGHT MANUAL WHEELCHAIRS

The V200 D is a light weight manually propelled wheelchair for temporary or permanent use with optimised settings.



C59
champagne



C60
brown

- 4 seat heights
- 2 seat depths
- 9 seat widths
- 3 armrest heights
- 3 positions for the armpads
- Mobile phone pouch under the seat
- A range of available accessories

FEATURES

A multitude of adjustments and a wide variety of options



Angelwings backrest, providing comfort for specific needs for your back



- Pushhandles adjustable in height
- String bag


In 3 heights and depths adjustable armpads

Foldable and detachable footrests

Rear wheels with Quickrelease system

Angle adjustable footplates

Wide range of accessories and options.



590	320	390	220-240	440-530	400	345-510	440-460	1010	16,5	130
620	320	420	220-240	440-530	400	345-510	440-460	1010	16,5	130
640	320	440	220-240	440-530	400	345-510	440-460	1010	16,5	130
660	320	460	220-240	440-530	400	345-510	440-460	1010	16,5	130
680	320	480	220-240	440-530	400	345-510	440-460	1010	16,5	130
700	320	500	220-240	440-530	400	345-510	440-460	1010	16,5	130



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.