

Sagitta

Lightweight active wheelchair

The Sagitta is a modern, active wheelchair that can be fully personalised. Armrests, leg supports and adjustable dimensions are included.




VERMEIREN
we care for you

Lightweight

The maximum weight is only 9.5 kg due to the solid frame made of light aluminium. Available in a wide range of colours, spoke protectors with colourful prints and a wide range of front and rear wheels.

Colours



Features



Backrest tilt adjustment
without tools



Available in SI-version
with reclining seat



Adjustable wheel camber
from 0° to 8°



Standard handbrakes

Technical data

Total width	550 - 780 mm
Total length	755 - 855 mm
Total height	720 mm
Folded width	–
Seat width	360 - 460 mm
Armrest height	90 mm

Seat height	470 - 520 mm
Backrest height	230 - 430 mm
Seat depth	330 - 380 - 420 mm
Weight	11,7 kg
Maximum user weight	100 kg
Additional information	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Sagitta** complies with the highest norms, laid down by certified test institutes in most European countries.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer