

# 268 - 269

## WALKING FRAMES

### Model 268 foldable with table

Tubular steel frame painted in blue.  
Padded arm support adjustable in height.  
4 wheels of 75 mm of which 2 with brake.  
Can be dismantled for easy transport and storage.



**C27**  
blue



### Model 269 Foldable with lower arm support

Tubular steel frame painted in blue.  
Lower arm support adjustable in height.  
4 wheels of 75 mm of which 2 with brake.  
Can be dismantled for easy transport and storage.

### Options for walking frames

- Table (type 268)
- Lower arm supports (type 269)
- Seat



268  
269

830-1130  
830-1130

760  
760

520  
520

580  
580

15,5  
12,5

120  
120

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

