



Material			
No.	Plan	Article	Qty
1	Y017225	1743754	4
2*	X001308	G051799	2
3*	X001309	G051899	2
4*	X002265	1500054	2
5**	X002437	1992699	2
6	X003635	1633399	8
7**	X003662	1998693	2
8**	X005720	1905243	2
9**	M16	1906194	2
10**	---	1991359	2
11**	M16 DIN17222	1911653	2
12	M6x45 DIN912	2994058	8
13	M6 DIN986	2994059	8
14	M6 DIN125	2994414	8
15*	XC00414	3991378	2
16*	M12x95 DIN931	3992074	2
17	E00052594	1919567 (L) 1919568 (R)	1 1

Quantities are for one wheelchair

* For fixed wheel

** For quick release wheel

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

V500 LIGHT T30			
Scale: 1:5	Drafted	Verified	Confirmed
Date	14/02/2023	22/03/2023	22/03/2023
Initials	SD	SD	PJ

Plan No.
E00052553/B_1.1

