



Turios Battery Harness Set 1 is used for:

1. Turios without lights and without electrical seat and backrest angle adjustment
2. Turios with electrical seat and backrest angle adjustment

Plan No.			
E00008074/A			
Scale:	1:1	Drafted	Verified
Date	04/11/2021	04/11/2021	04/11/2021
Initials	KS	MVG	TK



TURIOS BATTERY HARNESS SET 1 Sheet 1/2



No.	Plan	Article	Quantity
11	---	1994573	1
12	---	1901755	2
13	---	H005006	0,33

11 ← 550mm

13 ← 330mm

TURIOS BATTERY HARNESS SET 1 Sheet 2/2

Scale:	1:1	Drafted	Verified	Confirmed
Date	04/11/2021	04/11/2021	04/11/2021	04/11/2021
Initials	KS	MVG	TK	

Plan No.
E00008074/A